

CLEAR TRAINING CARTRIDGES

PHYSICIAN TRAINING OUTSIDE THE CATHLAB




Biomodex® Clear Training Cartridges are 3D printed using a transparent material that allows for training outside of the cathlab. Clear cartridges are inserted into Biomodex® EVIAS™ stations to simulate all ischemic stroke and unruptured aneurysm procedures.

BENEFITS



Anywhere/Anytime Flexibility



Quickly Swap Cartridges
for Different Anatomies



Varying Complexity and
Tortuosity for Iterative Training

WARNING: This product is not a medical device and is not intended to diagnose, treat, cure, or prevent any health condition. Biomodex® products have been designed for training and their safety and effectiveness as medical devices have not been demonstrated.

CLEAR TRAINING CARTRIDGES



TRAIN FOR THE MOST COMMON NEUROVASCULAR PROCEDURES INCLUDING:



Clear Cartridges - Standard Frames

Ischemic Stroke:
Mechanical thrombectomy
Aspiration thrombectomy

Unruptured Aneurysm:
Flow diversion
Intrasaccular device
Coil embolization
Balloon-assisted coiling
Stent placement

KEY FEATURES



Clear Cartridges - Standard & Tank Frames

Train for all unruptured aneurysm and ischemic stroke procedures
Transparent material designed for use without fluoroscopy
Compatible with all EVIAS™ stations
Standard frames for improved lateral visualization and convenience
Water-bath tank frames for improved transparency of the anatomy

COMPATIBILITY



EVIAS™ Station



EVIAS™ Plus Station



Thrombotech™ Synthetic Clot

To order

Call (813) 468-8086 or visit orders.biomodex.com

For more information

email contact@biomodex.com or visit biomodex.com